Wolfram Droh GmbH

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REV GA/2024-01

Yuyao Shengchang Medical Instrument Factory
Simen-Town
315400 Yuyao City, Zhejiang
China

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GB Contents

1	General	Page	12
2	Safety and Responsibilities	Page	12
3	Modification of the product	Page	13
4	Before using the product	Page	13
5	Using the product	Page	14
6	Product Maintenance / Reprocessing	Page	14
7	Accessories	Page	14
8	Replacement parts and repairs	Page	14
9	Disposal	Page	15
10	Technical Specifications	Page	15
11	Explanation of symbols	Page	16

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1 General

- The product may be used by trained personnel as well as by lay persons.
- The product is intended to be used to supply, inhale and humidify medical oxygen within homecare oxygen therapy as well as in medical settings.
- An oxygen mask with reservoir bag can be used to achieve the required high oxygen concentration.
The instructions of the physician are to be observed.
- An oxygen catheter is suitable for short-term use and up to a maximum of 4 lpm.
- The following instructions must be followed by untrained users.
- These instructions reflect the state of the art at the time of publication.
- This product complies with the provisions of the EU Directive 93/42/EEC (Medical Device Directive) and Medical Device Regulation (EU) 2017/745.
- The instructions for use comply with DIN EN 82079-1.

2 Safety and Responsibilities

IMPORTANT! READ CAREFULLY BEFORE USE
KEEP FOR FUTURE REFERENCE

- Please observe the following instructions and care instructions before use.
- If in doubt about the use and connection of the product, please contact customer service.
- If additional oxygen is administered care should be taken to ensure that smoking and other naked flames are forbidden in the vicinity of oxygen-carrying equipment and fittings.
- Do not smoke!
- The product draws air from the immediate surrounding and must therefore not be used in a toxic atmosphere.
- Keep the product free of oil and grease, otherwise there is a risk of explosion. This includes the exposure to cleaning alcohols, hand cream and adhesives – in contact with highly compressed oxygen explosive reactions may occur.
- The product contains di(2-ethylhexyl) phthalate (DEHP) and should therefore only be used temporarily in pregnant women, nursing mothers, babies and young children.
- A maximum length of 18 m (tubing and inhalation unit) should not be exceeded. Otherwise, the pressure and the oxygen flow cannot be guaranteed.
- The product must not be kinked, stretched or bent.
- Use a swivel connector to avoid twisting the tube.
- **Single use products may not be reused!**
They must be disposed of immediately after use in accordance with chapter 9.

The reuse of single use products represents a potential risk for the patient or the user. It can result in contamination and/or reduced functionality.

- Contamination and/or limited functionality of the product can lead to injury, illness or death of the patient.

3 Modification of the product

- Apart from the described maintenance instructions, the product may not be modified or repaired by the operator!
- The manufacturer's responsibility is limited to the dimensional and physiological suitability of the product.
- Improper use, modifications to the product or negligence will void any warranty. This also applies to unauthorised repairs. Normal signs of wear and tear on the device are also excluded.

4 Before using the product

- Do not expose the product to direct sunlight or other light sources. Store and transport the product in a dry and dust-free place in its original packaging at -10°C to $+50^{\circ}\text{C}$.
- Remove all delivered parts from the packaging. If necessary, remove the respective plastic covers and dispose of the packaging according to the current national regulations.
- Before each use verify that parts or components do not show any signs of hardening, damage, breakage or cracks. The product must also not be contaminated or soiled. If any of the conditions described above are present, the device must not be put into operation and must be disposed of in accordance with chapter 9.
- In case of doubt perform a comparison check with a brand-new product.
→ **Damaged products must not be used!**
After consulting customer service, the product must be disposed of according to chapter 9.
- The product contains PVC.
- Please note that the product is supplied in a non-sterile condition.
- The hygienic preparation of the product lies within the clinical responsibility of the user! The plastic product can be disinfected with a suitable cleaning agent (e.g. alcohol) before use.
Sterilization is not permitted!
- To attach and detach the tubing, the tubing must be held at the connector to prevent cracks from forming between the connector and the tubing.
- **Mounting the oxygen catheter:**
Connect the fittings to the corresponding adapters and/or connectors of the oxygen tubing. Insert the end of the catheter with the securely attached dressing into a nostril.
For total lengths over 15 m, a water trap should be used to collect any condensation.

- **Connecting the nasal cannula:**

Connect or screw the fittings to the appropriate connectors or directly to the oxygen source.

Insert the nasal prongs at the other end of the tubing directly into the nostrils. Place the ends of the tubing over the ears and secure these under the chin using the bolo.

This should be pulled up to the chin for security.

When using a nasal cannula with earpieces, insert the nose prongs into the nostrils and place the earpieces behind the ears.

- **Connecting the O₂ mask:**

Connect the fittings to the appropriate connectors or directly to the oxygen source.

Place the mask over mouth and nose onto the face. Pull the elastic over the head and under the ears at the nape of the neck. For a tight fit of the mask, tighten the elastics and adjust the metal nose piece to the patient's face.

- When completely assembled, check that all parts are connected tightly.
- The product is now ready for use.

5 Using the product

- The product must not come into contact with sharp objects or rub against an edge during use.
- After successful assembly, the product can be attached to the oxygen source. Please observe the requirements of the other connection devices and products as well as the instructions of the physician.
- The flow rate for inhalation systems is max. 6 lpm for the standard version.
- In addition, note the inspiratory oxygen concentration in relation to the flow under chapter 10 and set the flow as specified by the physician in order to obtain an optimum supply.

6 Product maintenance/Reprocessing

- Single use products may not be reused!
- The duration of use is limited to a maximum of 7 days.
In case of contamination or defects, replace the product immediately.
- The product must not be used on more than one patient.
- The nasal prongs/ the oxygen mask and the headpiece of the inhalation unit should be cleaned daily with a clean, slightly damp cloth.
- After termination of the application, the product must be disposed of according to chapter 9.

7 Accessories

- The product can be connected to other tubes via a connector.
- The user is responsible for the correct combination with other medical products.

8 Replacement parts and repairs

- In case of functional or material problems, contact customer service immediately.

9 Disposal

- The product must be disposed of after use.
- Please note the expiry date if applicable.
The product must not be used after the expiry date!
- Safe disposal is carried out in accordance with current national regulations.

10 Technical Specifications

Oxygen concentration in relation to the standard flow rate: 1 LPM = 4%		
Article	Oxygen flow (lpm)	Inspiratory oxygen concentration - FiO ₂
Nasal cannulas	1	21% - 24%
	2	24% - 28%
	3	28% - 34%
	4	34% - 38%
	5	38% - 42%
	6	42% - 46%
Oxygen mask	5 - 10	40% - 60%
Oxygen mask with reservoir, , without rebreathing	5 - 7	40% - 90%

Patient group	Age	Weight
Babies	up to 6 months	> 2,5 kg
Infants	3 - 24 months	> 4,0 kg
Children	1 - 6 years	> 12,5 kg
Small adults		30 - 70 kg
Adults		> 70 kg

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