DAR™

Guedel Airway

en

Guedel Airway

Intended Use/Purpose/Indications: The Guedel Airway is an oropharyngeal airway which is intended to be inserted into a patient's pharynx through the mouth to provide a patent airway. It is indicated to help restore airway patency.

Contraindications: Do not use on conscious patients or patients with a cough, gag reflex, or patients with a foreign body obstruction in the airway.

Warnings/Cautions: Do not clean and/or sterilize, single use only. Attempts to clean or sterilize result in bio-incompatibility, infection, or product failure risks. Device and its packaging are not made with natural rubber latex or DEHP. Only for use by trained clinicians. Inspect before use for damage or dirt. If damaged or the sterile barrier has been compromised do not use and contact Covidien. Damaged devices may result in an aspiration hazard. Replace the device if obstruction occurs from secretions or other material. Use duration should not exceed twenty-four (24) hours. Dispose of product in accordance with current medical standards and applicable national regulations for biologically hazardous waste.

Adverse Events: Adverse events associated with the procedural placement and use of this device in descending order of severity and frequency are: Tissue damage/tissue trauma, blunt trauma, respiratory distress, aspiration, infection, respiratory failure, delay of treatment, and anaphylaxis. Any serious incident that occurs in relation to the device or its accessories should be reported to Covidien and to the national competent authority.

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	Sterilized using ethylene oxide
	Non-sterile
	Single use
	Do not use if package is damaged
	Consult instructions for use
	Caution
	Keep away from sunlight
	Keep dry
	Fragile
	Authorized representative in the European Community
	Catalogue number
	Manufacturer
	Use-by date
	Lot number
	Date of manufacture

Part No. PT00113490 Rev A 2020-01

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG.

TM brands are trademarks of their respective owners. Other brands are trademarks of a Covidien company.

© 2019 Covidien.
Covidien Inc, 15 Hampshire Street, Mansfield, MA 02048 USA.
EC REP Medtronic B.V., Earl Bakkenstraat 10, 6422 PJ Heerlen, The Netherlands.

www.covidien.com

Rx
ONLY

CE
0123