
McGRATH™

Rechargeable Battery System

en	Instructions For Use
fr	Mode d'emploi
de	Gebrauchsanweisung
nl	Gebruiksaanwijzing
it	Istruzioni per l'uso
es	Instrucciones de uso
sv	Bruksanvisning
da	Brugsanvisning
no	Bruksanvisning
fi	Käyttöohjeet
pt-BR	Instruções de uso
pt-PT	Instruções de utilização
ru	Инструкция по применению
pl	Instrukcja użytkowania
cs	Návod k použití
sk	Pokyny na používanie
sl	Navodila za uporabo
hu	Használati útmutató
el	Οδηγίες χρήσης
tr	Kullanım Talimatları
bg	Инструкции за употреба
ro	Instrucţiuni de utilizare
et	Kasutusjuhend
lt	Naudojimo instrukcija
lv	Lietošanas pamācība
hr	Upute za upotrebu
sr	Uputstvo za upotrebu
id	Petunjuk Penggunaan

REF 341-100-000**REF** 341-200-000

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IMPORTANT! Recharge the battery to full capacity prior to first use.

Description: The McGRATH™ rechargeable battery system is a power source for the McGRATH™ MAC video laryngoscope including the rechargeable battery and multi-battery charger. The rechargeable battery powers the video laryngoscope handle during use. The McGRATH™ multi-battery charger can wirelessly recharge five batteries simultaneously.

Intended use: The McGRATH™ rechargeable battery is a power source to McGRATH™ MAC video laryngoscope which is intended to be used by trained and licensed individuals to view the vocal cords during medical procedures, typically laryngoscopy.

Intended purpose: The McGRATH™ rechargeable battery powers the McGRATH™ MAC video laryngoscope during visualization of the oropharynx and larynx.

Intended user: The McGRATH™ rechargeable battery is intended for use by trained and licensed individuals during anesthesia, intensive care and emergency care. Intended users include:

- Nurse anesthetists (CRNA) or equivalent clinical role
- Respiratory therapists (RT) or equivalent clinical role
- Anesthetists or equivalent clinical role
- Emergency department physicians or equivalent clinical role
- Physicians, critical care

Target patient population: The McGRATH™ rechargeable battery can be used on pediatric and adult patients that require a view of the vocal cords during medical procedures.

Indications: General laryngoscopy and to aid in securing the airway.

Contraindications: None known.

Catalogue number	Product descriptions
341-100-000	McGRATH™ rechargeable battery
341-200-000	McGRATH™ multi-battery charger
341-400-000 (Accessory)	McGRATH™ multi-battery charger cover
341-500-000 (Accessory)	McGRATH™ multi-battery charger mounting bracket
341-600-000 (Accessory)	McGRATH™ multi-battery charger power supply and adapters

For ordering information, please contact your local representative.

Warnings

- Only personnel trained in and licensed to perform intubation with a laryngoscope may use this device.
- Do not use accessories other than those specified or provided by the manufacturer.
- Do not open any sealed part of the rechargeable battery or multi-battery charger. Doing so will impair the performance of the device and void the warranty, and it may impact patient safety.
- Do not use the rechargeable battery or multi-battery charger if any components become loose, physically damaged, or fit poorly.
- Do not use the rechargeable battery or multi-battery charger in the presence of flammable mixtures.
- Follow local regulations for disposal of the rechargeable battery and multi-battery charger when it reaches the end of life.
- Do not use the rechargeable battery or multi-battery charger in a magnetic resonance (MR) environment.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Portable radio-frequency communications equipment (including peripherals such as antenna cables and external antennas) should not be used within 30 cm (12 in) of the multi-battery charger as this can degrade performance.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer, could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- The multi-battery charger meets the requirements for radio-frequency exposure when operated at a minimum distance of 20 cm from the user or nearby persons.
- Failure to adhere to stated warnings may lead to a fire or explosion with the rechargeable battery, potentially causing a severe burn injury.
- Ensure the rechargeable battery is completely dry before use.
- Remove the rechargeable battery from the charger prior to cleaning or disinfection.
- Do not autoclave the rechargeable battery.
- Do not reprocess the rechargeable battery in an ultrasonic cleaner.
- Do not dispose of the rechargeable battery in a fire. Do not crush, short circuit, or expose the battery to temperatures above 60°C (140°F).
- Do not clean or disinfect these devices using any method other than what is listed in the "Cleaning and Disinfection" section of this document. Doing so may damage these devices and lead to delay of treatment.
- Do not allow the rechargeable battery to be mechanically shocked or become deformed by external forces.
- Do not use the rechargeable battery after the expiry date.

Cautions

- Federal (USA) law restricts the sale of this rechargeable battery to purchase by or on the order of a licensed practitioner.
- Do not put any other type of rechargeable battery into the multi-battery charger.
- Do not put any foreign objects in the electro-magnetic field of the multi-battery charger, as they might become hot.
- Disconnect the power supply from the multi-battery charger before cleaning.
- Do not store the rechargeable battery long term in its fully charged state. Doing so may result in partial loss of battery capacity and permanent degeneration in next charging cycle.
- Install the multi-battery charger at a position where the AC/DC adapter can easily be disconnected from the AC power supply if an electrical problem occurs, otherwise the multi-battery charger may be damaged due to electrical problems.
- The rechargeable battery is designed for use with the McGRATH™ MAC video laryngoscope only.
- The maximum storage time of the rechargeable battery is 3 months. At least once every 3 months the battery needs to be discharged to 0 minutes, and then fully recharged to 120 minutes. Failing to do so will decrease or even end the battery life.
- The solid purple LED on the rechargeable battery means the battery is in error mode, or it has reached the end of life and should be disposed of in compliance with local regulations.
- The VHB tape is designed for mounting on the wall flat with the unit affixed to the bracket. Ensure no bubble, wrinkle or uneven surfaces are on the tape. Follow the assembly sequence when mounting the bracket on charger.
- The bracket screw is for the table stand application only. Follow the multi-battery charger mounting bracket instructions for proper assembly.

- The charger can support at most five rechargeable batteries simultaneously (5-2). If over five batteries are placed into the charger, the batteries may not be charged properly.
- Be careful when mounting the charger on the wall and removing it. Dropping the charger could result in charger malfunction.

Adverse events: The following adverse events have been reported in association with the use of the rechargeable battery during laryngoscopy and intubation. Listed in order of severity the adverse events are: exposure to body fluids, unspecified tissue injury (including blunt trauma), respiratory failure (including hypercarbia), hypoxia, unspecified infection, burn(s), electrical shock, delay of treatment (prolonged procedure/extended hospitalization), chemical exposure (including toxicity), unintended radiation exposure. If a serious incident related to the device occurs, immediately report the incident to Covidien and the applicable competent authority.

Specifications

McGRATH™ rechargeable battery	
Size:	85 mm x 32 mm x 19 mm
Weight:	25 g
Battery type:	Rechargeable Li-ion Battery
Capacity:	520 ±20 mAh, 1.97 Wh
Protection:	IP33
Power:	2 W
Display:	Color LED indicator
Nominal voltage:	3.8 VDC
Output voltage:	3~4.35 VDC
Electromagnetic emissions:	CISPR 11 Class B, Group 1
Electromagnetic immunity:	IEC 61000-4-2, IEC 61000-4-3, IEC 61000-4-8
Useful life:	4 years or a maximum of 10000 minutes operating time
McGRATH™ multi-battery charger	
Size:	245 mm x 133 mm x 46 mm
Weight:	845 g
Power supply:	Proprietary McGRATH™ power supply
Protection:	IP22
Charging output power:	9 W (Max)
Display:	Color LED indicator
Electromagnetic emissions:	CISPR 11 Class A, Group 2 IEC 61000-3-2 Class A, IEC 61000-3-3
Electromagnetic immunity:	IEC 61000-4-2, IEC 61000-4-3, IEC 61000-4-4, IEC 61000-4-5, IEC 61000-4-6, IEC 61000-4-8, IEC 61000-4-11, RFID tag interference test (1 mm to 25 mm)
Useful life:	4 years
VHB double sided tape	
Useful life:	2 years
Power supply	
Input voltage and frequency:	100~240 VAC, 50 Hz / 60 Hz
Output voltage and power:	15 VDC, 18 W
Interchangeable AC adapter plugs:	 Note: Choose the appropriate plug per country.

Environmental Conditions

Operating and transient conditions	
Device temperature	10 °C (50 °F) to 40 °C (104 °F)
Atmospheric pressure	620 hPa to 1060 hPa
Relative humidity	10% to 95% (non-condensing)
Storage and transportation	
Device temperature	-10 °C (14 °F) to 40 °C (104 °F)
Atmospheric pressure	620 hPa to 1060 hPa
Relative humidity	10% to 95% (non-condensing)

Unpacking, inspection: The McGRATH™ rechargeable battery system is comprised of 5 rechargeable batteries at most, 1 multi-battery charger, 1 charger cover, 1 bracket, and 1 power supply. Unpack and check for any visible signs of damage. In the event of damage, do not use the device and inform your local representative. See Figure 2 in the Quick Start Guide.

Operation: Always check the battery power remaining in the McGRATH™ MAC video laryngoscope before starting a procedure. The power is given in minutes on the handle display as shown in Figure 9 of the Quick Start Guide.

Battery: The rechargeable battery can power McGRATH™ MAC video laryngoscope 301-000-000 handle. Push the battery unit into the cavity of the handle as shown in Figure 5 of the Quick Start Guide. When the battery is secured in the McGRATH™ MAC video laryngoscope, the rechargeable battery can be turned on or off by pressing the power button.

The battery has an LED indicator that shines through the battery housing and provides a visual indication of how much power remains in the battery. See the "LED indicator status" section for more information. The battery's LED indicator only functions when not installed in a video laryngoscope.

When the battery is placed in a McGrath™ MAC video laryngoscope handle and the handle is turned on, the remaining operating time of the battery unit in minutes is displayed in the bottom right corner of the video laryngoscope screen. See Figure 6 in the Quick Start Guide. This counts down as battery power is consumed. The operating time on a full charge is 120 minutes.

When the battery counter reaches 5 minutes on the 301-000-000 video laryngoscope screen, the battery icon will flash. If the counter is allowed to run down to zero, the battery icon will continue to flash and the battery supports at least 10 minutes of operating time before the 301-000-000 video laryngoscope is shut off automatically. Or, the user can press the battery button to switch off the video laryngoscope and pull on the tab to remove the battery.

Charger: The charger is powered by an AC/DC power supply. Prior to each use, check the integrity of the power supply. Peel off the instructional label on the charger cover before the first time use. The charger can be mounted on a wall, placed on a table, or laid on a flat surface. See the Quick Start Guide for assembly instructions for each of these orientations. See bracket Instructions For Use for the usage of VHB tape and bracket screw during the wall-mount and table-stand assembly.

- Wall-mount assembly (Figure 6 in the Quick Start Guide):

- (6-1) Use a 70% isopropyl alcohol (IPA) wipe to clean the bracket surface and wait for it to dry. (6-2) Pay attention to the bracket orientation as shown in the figure.
- Peel off the bottom liner of VHB tape and affix the tape on the designated area of bracket. (6-3) Press the tape firmly for 15 seconds. Ensure the VHB tape is applied evenly with no bubbles or wrinkles.
- (6-4) Use 70% isopropyl alcohol (IPA) wipe to clean the wall surface (coated sheet metal or painted drywall) and wait for it to dry. (6-5) Peel off the front liner of VHB tape. (6-6) Turn the bracket around and affix it on the wall by pressing the bracket firmly for 15 seconds. (6-7) Wait at least 24 hours to ensure maximum adhesive strength.
- Mount the charger on the bracket using the receiver on the back of the charger. (6-8) You will hear a snap if correctly assembled.
- (6-9) Insert the charger cover in the guides on the sides of the charger and slide it down.
- (6-10) Insert 1 to 5 batteries at most.

- Table-stand assembly (Figure 7 in the Quick Start Guide):


- (7-1) Turn bracket so that "McGrath" is facing down then mount the bracket to the slots on the bottom of charger. Rotate clockwise to tighten the bracket screw underneath.
- (7-2) Insert the charger cover in the guides on the sides of the charger and slide it down.
- (7-3) Insert 1 to 5 batteries at most.

- Lay-flat assembly (Figure 8 in the Quick Start Guide):

- Mount the bracket onto the back side of the charger. (8-1) Pay attention to the bracket orientation as shown in the figure. (8-2) You will hear a snap if correctly assembled. (8-3) Then check to make sure the bracket is mounted firmly without any movement.
- (8-4) Insert the charger cover in the guides on the sides of the charger and slide it down.
- (8-5) Insert 1 to 5 batteries at most.

Battery charging: Power the multi-battery charger prior to charging batteries. To charge the rechargeable battery, insert it into the multi-battery charger with the removal tab facing up. The charger could support at most 5 rechargeable batteries charging. Once in the charger, the rechargeable battery can be fully charged in 6-8 hours. **Note:** (5-1) Ensure the battery was cleaned and disinfected prior to charging.

Battery maintenance: If the rechargeable battery remaining in the video laryngoscope handle is not in use for more than 1 month, remove the rechargeable battery from the video laryngoscope handle before storage. If the battery has not been used for more than 3 months, recharge the battery before placing it in the McGrath™ MAC video laryngoscope handle. The maximum storage time of the rechargeable battery is 3 months. At least once every 3 months in long-term storage, the battery needs to be discharged to 0 minutes, and then fully recharged to 120 minutes.

LED indicator status: The indicators reflect the status of the rechargeable battery and multi-battery charger. The rechargeable battery LED indicator has 4 colors to signal the corresponding status of the battery as shown in Figure 3 of the Quick Start Guide.

- Red LED: The battery has less than 25% power remaining. It can support the McGrath MAC video laryngoscope for less than 30 minutes.
- Amber LED: The battery has more than 25% power remaining. It can support the McGrath MAC video laryngoscope for more than 30 minutes.
- Green LED: The battery has more than 95% power remaining. It can support the McGrath MAC video laryngoscope about up to 120 minutes.
- Purple LED: The battery is in error mode, or it has reached the end of life. For both cases, the battery should be disposed of in compliance with local regulations.


The multi-battery charger LED indicator has 3 colors that provide a visual indication of the status of the charger as shown in Figure 4 of the Quick Start Guide.

- No LED: The charger is not powered.
- White LED: The charger is powered and is working normally.
- Red LED: The charger is not working due to an error. **Note:** Once the error is resolved, the charger returns to normal mode (white LED) without user intervention.

Cleaning and disinfection: The McGrath™ multi-battery charger and rechargeable battery are supplied as non-sterile. After unpacking, clean and decontaminate the charger and battery before use. The battery should be separated from the video laryngoscope after each patient use, and cleaned and disinfected per the instructions for use and in accordance with local regulations.

The following instructions have been validated by Covidien as adequate for preparing the device for reuse.

Approved cleaning and disinfection methods:

Product	Cleaning	Disinfection
McGRATH™ multi-battery charger 	Use 70% isopropyl alcohol (IPA) wipe, or quaternary ammonium + isopropyl alcohol (IPA) wipe (for example, Super Sani-Cloth or equivalent).	Not required.
McGRATH™ rechargeable battery 	Use 70% isopropyl alcohol (IPA) wipe, or quaternary ammonium + isopropyl alcohol (IPA) wipe (for example, Super Sani-Cloth or equivalent).	1. Clean. 2. Using quaternary ammonium + isopropyl alcohol (IPA) wipe (for example, Super Sani-Cloth or equivalent) for disinfection.

Dispose of wipes in accordance with local regulations and hospital policy.

Multi-battery charger cleaning

To clean the charger:

- (6-11) First, remove the batteries and charger cover, then remove the power supply from the charger. **Note:** If the charger is mounted on the wall, after emptying batteries, remove the charger cover. (6-12) Then detach the charger from the bracket on the wall by pressing the bracket latch and lifting the charger from the bracket.
- Use cleaning material such as a 70% isopropyl alcohol (IPA) wipe, or quaternary ammonium + isopropyl alcohol (IPA) wipe (for example, Super Sani-Cloth or equivalent) to clean the charger cover and charger housing as described in Step 2a through Step 2b. Ensure a minimum contact time of 1 minute.
 - Charger cover: Treat all surfaces with the cleaning material to obtain thorough coverage.
 - Charger housing: Clean all cabinet surfaces.
- Visually inspect the charger for visible soil. If soil is still visible, use a new wipe to repeat the cleaning process, beginning with Step 2.

Rechargeable battery cleaning and disinfection

To clean the battery:

Use cleaning material such as 70% isopropyl alcohol (IPA) wipe, or quaternary ammonium + isopropyl alcohol (IPA) wipe (for example, Super Sani-Cloth or equivalent). Clean the battery from top to bottom, systematically working through the following 5 steps. Apply the cleaning material to all surfaces of the battery. Ensure a minimum contact time of 1 minute.

- Power contact: Apply cleaning material to each hole and corner of the metal battery contacts. Scrub gently, taking care not to deform or damage the metal terminals.
- Retaining clips: Run the cleaning material along and into the gaps around the battery retaining clips.
- Lower battery slot: Scrub every hole and corner surface of the rechargeable battery.
- Battery housing: Wipe all surfaces of the rechargeable battery housing with the cleaning material.
- Removal tab: Run the cleaning material along all tab surfaces.

If soil is still visible, use a new wipe to repeat the cleaning process. Ensure that the cleaning material contacts all surfaces of the device even if no visible soiling is present.

To disinfect the battery:

First, clean the battery per the instructions above. Then use a new quaternary ammonium + isopropyl alcohol (IPA) wipe (for example, Super Sani-Cloth or equivalent) to disinfect the battery, systematically working through the same steps for cleaning the battery.

Storage: After cleaning and disinfection, air dry or wipe the multi-battery charger and rechargeable battery with a lint-free cloth. Dry the charger and batteries completely before storage. Put the decontaminated charger and batteries in a storage container immediately after drying and store according to local guidelines. The device and batteries cannot be stored in the original packaging for reuse.

Troubleshooting

McGRATH™ multi-battery charger:

Scenario	Possible causes	How to resolve
Charger LED red	Interfering object (for example, metal object) in the charging area	Remove interfering object (for example, metal object). Charger will recover automatically after 10 seconds.
	Charger overheats	Wait until charger temperature drops to below 60 °C
	Abnormal working voltage	Check power supply connection
Charger LED off	Wrong power supply connected or power cord malfunction	Replace the power supply and restart the system

McGRATH™ rechargeable battery:

Scenario	Possible causes	How to resolve
Battery LED off in the charger	Battery not in the right position for charging	Ensure the battery tab is facing up and the power button is facing outwards
	Charger in error mode	Ensure the charger is powered on and working normally
Battery LED solid purple	Battery reached end of life or permanent failure	Dispose of battery according to local regulations and purchase a new battery
Battery LED alternating purple	600 minutes of useful life remaining	Consider replacing the battery
No image shown on video laryngoscope screen	Battery power depleted	Recharge the battery
	Physical damage (for example, power contact broken, battery button malfunction)	Replace with new battery
Battery failed to reactivate by pressing on button or recharging	Battery capacity exhausted or malfunction	If the battery cannot be activated due to battery damage or lost efficacy, replace with new battery.

If none of the troubleshooting measures achieve positive results, return the device or devices to your Covidien agent or distributor for a replacement.

Electromagnetic Interference (EMI)

The McGRATH™ rechargeable battery system has been tested to product specific standard IEC 60601-1-2:2014+A1:2020 (Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests).

The following tables contain the corresponding manufacturer's declaration for the electromagnetic emissions and electromagnetic immunity of the device.

Electromagnetic Compatibility (EMC)

The McGRATH™ rechargeable battery system is intended for use in professional healthcare facility environments.

The McGRATH™ rechargeable battery system complies with international standards for radiated and conducted emissions, susceptibility to radiated emissions, electrostatic discharge, and power disturbances according to the applicable standards.

Cable	Maximum cable length, shielded/unshielded	Number	Cable classification
AC power line	1.8 m, unshielded	1 set	DC power

The McGRATH™ rechargeable battery system requires special precautions during installation and operation for electromagnetic compatibility. In particular, the use of nearby mobile or portable communications equipment may influence the performance of the McGRATH™ rechargeable battery system.

Electromagnetic emissions

Guidance and Manufacturer's Declaration - Electromagnetic Emissions	
Emissions test	Compliance
RF emissions CISPR 11	Group 1 (Rechargeable battery)
	Group 2 (Multi-battery charger)
RF emissions CISPR 11	Class B (Rechargeable battery)
	Class A (Multi-battery charger)
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies

Electromagnetic immunity

Guidance and Manufacturer's Declaration - Electromagnetic Immunity		
Immunity test	IEC/EN 60601-1-2 test level	Compliance level
	Electromagnetic environment guidance	
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air
Electric fast transient/burst IEC 61000-4-4	± 2kV (100 kHz repetition rate) AC mains ± 1kV (100 kHz repetition rate) input/output >3 m	± 2kV (100 kHz repetition rate) AC mains ± 1kV (100 kHz repetition rate) input/output >3 m
	Mains power quality should be that of a typical commercial or hospital environment.	
Surge IEC 61000-4-5	± 1kV line-line, AC mains	± 1kV line-line, AC mains
	Mains power quality should be that of a typical commercial or hospital environment.	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	100% reduction for: • 0.5 cycle (at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°) • 1.0 cycle (at 0°) • 250/300 cycles (at 0°) 30% reduction for: • 25/30 cycles (at 0°)	100% reduction for: • 0.5 cycle (at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°) • 1.0 cycle (at 0°) • 250/300 cycles (at 0°) 30% reduction for: • 25/30 cycles (at 0°)
	Mains power quality should be that of a typical commercial or hospital environment.	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	50 and 60 Hz, 30 A/m, x-, y-, and z-axes	50 and 60 Hz, 30 A/m, x-, y-, and z-axes
	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.	
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz; 6 Vrms, ISM and amateur radio bands between 0.15 MHz and 80 MHz frequencies	3 Vrms 150 kHz to 80 MHz; 6 Vrms, ISM and amateur radio bands between 0.15 MHz and 80 MHz frequencies
	Except as indicated in the table "Test Specifications for Enclosure Port Immunity to RF Wireless Communications Equipment", portable and mobile RF communications equipment should be used no closer to any part of the McGRATH™ Rechargeable Battery System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ 150 kHz to 80 MHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey¹, should be less than the compliance level in each frequency range². Interference may occur in the vicinity of equipment marked with the following symbol: 	

Guidance and Manufacturer's Declaration - Electromagnetic Immunity		
Immunity test	IEC/EN 60601-1-2 test level	Compliance level
	Electromagnetic environment guidance	
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m 80 MHz to 2.7 GHz
	Except as indicated in the table "Test Specifications for Enclosure Port Immunity to RF Wireless Communications Equipment", portable and mobile RF communications equipment should be used no closer to any part of the battery system, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz and $d = 2.3 \sqrt{P}$ 800 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ¹ , should be less than the compliance level in each frequency range ² . Interference may occur in the vicinity of equipment marked with the following symbol: 	
RF wireless proximity fields IEC 61000-4-3	See table "Test Specifications for Enclosure Port Immunity to RF Wireless Communications Equipment". 0.3 m	
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.		
- Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the McGRATH™ rechargeable battery system is used exceeds the applicable RF compliance level above, the McGRATH™ rechargeable battery system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the McGRATH™ rechargeable battery system. - Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.		

Recommended separation distances between portable and mobile RF communications equipment and the McGRATH™ rechargeable battery system				
The McGRATH™ rechargeable battery system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the McGRATH™ rechargeable battery system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the McGRATH™ rechargeable battery system as recommended below, according to the maximum output power of the communications equipment except as indicated in the table "Test Specifications for Enclosure Port Immunity to RF Wireless Communications Equipment".				
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)			
	$d = 1.2 \sqrt{P}$ 150 kHz to 80 MHz out of ISM band and amateur band	$d = 0.59 \sqrt{P}$ 150 kHz to 80 MHz in ISM band and amateur band	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz	$d = 2.3 \sqrt{P}$ 800 MHz to 2.7 GHz
0.01	0.12	0.059	0.12	0.23
0.1	0.38	0.19	0.38	0.73
1	1.2	0.59	1.2	2.3
10	3.8	1.9	3.8	7.3
100	12	5.9	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.				
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.				
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.				

Proximity magnetic fields		
Test frequency	Modulation	Immunity test levels
		Professional healthcare facility environment
30 kHz	CW	8 A/m
134.2 kHz	Pulse modulation 2,1 kHz	65 A/m
13.56 MHz	Pulse modulation 50 kHz	7.5 A/m

Test Specifications for Enclosure Port Immunity to RF Wireless Communications Equipment. Frequencies and services of portable and mobile transmitters for which the recommended separation distance is 30 cm (12 inches)						
Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380 to 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430 to 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28

Test Specifications for Enclosure Port Immunity to RF Wireless Communications Equipment. Frequencies and services of portable and mobile transmitters for which the recommended separation distance is 30 cm (12 inches)

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
710	704 to 787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
745						
780						
810	800 to 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
870						
930						
1720	1700 to 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1845						
1970						
2450	2400 to 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5240	5100 to 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5500						
5785						

Warranty

The device is supplied with a manufacturer's warranty.

- Only products supplied by an approved Covidien agent or distributor are covered by the manufacturer's warranty.
- To be covered by the warranty, the product must be maintained in accordance with the procedures in this document.
- Before returning the device under a warranty claim, disinfect the device per the instructions in this document.

For the full terms of warranty, please contact your Covidien agent or distributor.

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COVIDIEN, COVIDIEN with logo, and the Covidien logo are trademarks of Covidien.
^{TM*} Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Covidien company.
© 2024 Covidien. All Rights Reserved.
Covidien Inc, 15 Hampshire Street, Mansfield, MA 02048 USA.
Medtronic B.V., Earl Bakkenstraat 10, 6422 PJ-Heerlen, The Netherlands.
Medtronic Limited, Building 9, Croxley Park Hatters Lane, Watford WD18 8WW, United Kingdom.