# LogiCal® Pressure Transducers

**en**

INSTRUCTIONS FOR USE

**fr**

MODE D'EMPLOI

**de**

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**it**

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**es**

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**pt**

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**nl**

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**sv**

BRUKSANVISNING

**da**

BRUGSANVISNING

**no**

BRUKSANVISNING

**fi**

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**el**

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# LogiCal® Pressure Transducer

## INSTRUCTIONS FOR USE

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

**NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.**

The Smiths Medical LogiCal® pressure transducer has been manufactured to meet or exceed all AAMI\* specifications for blood pressure transducers. The LogiCal® transducer is designed for use with Smiths Medical disposable monitoring kits and accessories. Smiths Medical disposable domes provide an isolation between fluid contacting the patient and the transducer diaphragm. Some Smiths Medical disposable monitoring kits are packaged with a flush device that is designed for use in pressure measurements which require continuous flow to maintain catheter patency (Fig. 3). The Smiths Medical flush device provides a controlled and constant fluid infusion at a rate determined by the patient's blood pressure and the fluid source pressure. With a total fluid source pressure equal to 300 mmHg, the flush device delivers a nominal flow rate at 3 ml/hour. (The neonatal flush, as indicated by an orange colour, delivers a nominal flow rate equal up to 30 ml/hour, the rate is determined by the setting on the syringe pump delivery system). Smiths Medical also offers disposable pressure monitoring without a flush device. This product may come in the form of a stand-alone dome or as an integral part of a complete monitoring kit used in cases where continuous flow is not a requirement.

\*Association for the Advancement of Medical Instrumentation

INTENDED USE

The Smiths Medical LogiCal® pressure transducer is a hemodynamic pressure monitoring system that is intended for blood pressure monitoring and blood sampling.

Contraindications: None known.

I. INTRODUCTION

1.

The transducer should be stored at ambient temperature. Storage conditions should never exceed a -25°C to +70°C range.

2.

Your LogiCal® transducer has been supplied ready to use. Additional mounting accessories are as follows:
- MX 260

MX 261P1

MX 262P1

MX 263P1

MX 264

MX 960XYP1

Pole clamp for mounting a transducer

Pole clamp for mounting a transducer or transducer brackets

Bracket for holding two transducers

Bracket for holding three transducers

Bracket for holding four transducers

Disposable dome
3.

The LogiCal® may be cleaned (see section V on Cleaning).

4.

The LogiCal® should be stored in a dry location. Transducers should be checked using a known pressure source prior to use or according to the specific monitor's instructions, and hospital protocol (see section IV).

5.

Connect the LogiCal® transducer to the monitor adaptor cable by mating the reusable cable connector with the transducer connector (i.e. plug), as shown in Fig. 2.
- PRECAUTIONS:
1.

Do not use any Smiths Medical pressure monitoring kits or transducers with non-isolated pressure monitors.
2.

Do not kink cables.

3.

Do not autoclave or use EtO to sterilize transducers, interface cables, clamps or brackets.

4.

Do not submerge the LogiCal® transducer.

5.

Allow a five (5) minute warm-up before calibration of system.

6.

Do not use a diaphragm that is cut or damaged. Use of a damaged diaphragm can result in inaccuracy of readings. This device must be thoroughly inspected for signs of damage prior to use.

7.

New split septum sampling site. Only use with needle free access devices; MX920NFC2, MX921NFC2, MX922NFC2, MX923NFC2.

8.

The disposable LogiCal® Pressure Monitoring kits should be changed every 96 hours or according to hospital protocol and/or CDC guidelines.
- WARNINGS:
- I.

1.

Do not use stopcocks or sampling site for the injection of fluids. Introducing pressure into the line (i.e. injecting fluids via syringe) may cause damage to the sampling site.

2.

This product contains DEHP – data derived from animal and epidemiological studies outline the following risks, but are not limited to: Adverse effects on fetal development, kidney and liver toxicity, negative impact on fertility. Special attention on the application of phthalate products should be given for the following patient populations: pediatrics, pregnant and nursing women. For further information on phthalates please refer to the Smiths Medical website: [www.smiths-medical.com/phthalates](http://www.smiths-medical.com/phthalates)

II.

PRIMING INSTRUCTIONS AND TIPS

A.

Priming with infusion bag:

PRIME WITH GRAVITY! - Applying pressure from pressure cuff may cause fluid to squirt and foam as it flows through the system. Thousands of microbubbles are formed and are difficult to remove. If fluid runs too slowly through the system, simply raise the supply bag.

Note: 1 cm of height is roughly 0.8 mmHg (every inch of height is roughly equal to 2 mmHg). A primed 100 cm line applies almost 80 mmHg of pressure if fully extended (a 4-foot line applies almost 100 mmHg of pressure if fully extended).

Note: The solution bag must be placed at a higher level than the transducer and pressure tubing in order to prime with gravity.

Step 1: Attach the dome to the transducer by sliding the dome down into the slots over the diaphragm until it clicks into place. Use this product dry-coupled only, no need to wet couple.

Step 2: Spike a bag of sterile solution with a gentle twisting motion of the drip chamber spike.

Step 3: Fill part of the drip chamber, by squeezing the sides, before starting to run fluid into the administration set.

Step 4: Open the roller clamp at the administration set and fill up to the flush device.

Step 5: Actuate the flush device gently (Fig. 3). Remember, air rises, so always allow the fluid to fill from below. Continue priming without pressure.

Step 6: Prime the side port of the zeroing stopcock and cap with a Luer lock plug.

CAUTION: Always turn the stopcock off to the transducer or off to the venting path before attaching the Luer Lock plug.

Step 7: Check the monitoring system for bubbles. Gently, tap areas of the kit where bubbles may be hidden from view.

CAUTION: Significant distortion of the pressure waveform or air embolism can result from air bubbles in the system.

Step 8: Pressurise the system with a Smiths Medical C-Fusor®, or Clear-Cuff® to 300 mmHg. Fast flush the system for 2-3 seconds, observe the device.
- 8a.

During flushing, fluid should flow freely from the end of the set and fluid flow will be visible in the drip chamber.

8b.

After flushing, fluid should not be visibly flowing from the end of the set and flow in the drip chamber should lie between the values listed in the table below:

Note:

The flow rate is affected by the number of lines being used.
- | Number of Lines in Set | Expected Flow Rate                                                                                   |
|------------------------|------------------------------------------------------------------------------------------------------|
| 1                      | Between 2 drops per minute and 7 drops per minute for the upper end of the performance specification |
| 2                      | 12 drops per minute will be equivalent to the upper end of the performance specification             |
| 3                      | 14 drops per minute are equivalent to the upper end of the performance specification                 |
| 4                      | 18 drops per minute are equivalent to the upper end of the performance specification                 |
- CAUTION: If fluid doesn't flow freely during flushing or if after flushing a higher flow rate is observed, the device is not functioning correctly and should not be used.

8c.

Check again for air bubbles which may have been dislodged during the fast flush.

Step 9:

The system is now ready to zero and calibrate.

B.

Priming with syringe:

Step 1: Attach a syringe filled with sterile saline solution to the system.

Step 2: Slowly depress the syringe plunger and flush with sterile saline solution until bubble free.

Step 3: After priming connect the transducer and desired tubing to the catheter.

Step 4: The system is now ready to zero and calibrate.

Note: Always position stopcocks at 90° for the "OFF" position. Do not place any stopcock handle at 45° to achieve an "OFF" position. Improper orientation of the stopcock handles may result in patient contamination, bleed back, or air embolism.

CAUTION: Avoid flushing air bubbles or clots, contained within the three-way venting stopcock or cannula, back into the patient by making sure that the monitoring line is completely filled with liquid, and by allowing a small amount of blood to flow back through the cannula before making the monitoring line connection.
- III. ZEROING
- Step 1:

Mount the transducer and disposable kit in the appropriate operating position with the stopcock zeroing port at the level of the right atrium, mid-axillary or 4th intercostal space.

Step 2:

Vent the transducer to atmospheric pressure by turning the zeroing stopcock off to the patient.

Step 3:

Adjust the monitor to read zero according to the monitor manufacturer's instructions.
- IV. CALIBRATION
- One of the following methods should be used to calibrate the monitor for sensitivity:

**Method A:** Apply a known pressure using a manometer or other pressure source, eg the check device 650-900 (use sterile technique and hospital protocol), and adjust gain or sensitivity according to the monitor manufacturer's instructions.

**Method B:** If supplied, press the "cal" button on the monitor or the 100 mmHg button on the transducer to check the calibration of the monitor. Adjust the gain or sensitivity according to the monitor manufacturer's instruction.
- Calibration will only be as good as the accuracy of the specified calibration requirement listed in the specifications section of these instructions.

V.

CLEANING

1.

Clean all blood and other foreign matter from the external surfaces of the transducer and cable using a slightly damp cloth and a mild detergent solution. Do not immerse. Do not use alcohol or solvents on the clamps.

2.

Rinse thoroughly.

WARNING! Measurement artifacts are possible with any pressure monitoring devices. Do not make changes to patient therapy based solely on the readings obtained from this device. Check that the patients' physical condition is consistent with any readings, failure to do so could result in inappropriate treatment.

CAUTION: Please tighten all connections before use!

DO NOT AUTOCLAVE OR ETO STERILISE THE LogiCal® TRANSDUCER.

LogiCal, Sure-Cal, C-Fusor, Clear-Cuff, Medex design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.

© 2019 Smiths Medical. All rights reserved.
- 2
- 3

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

| EN Specifications                                                                                                                            |  |                                                                           |
|----------------------------------------------------------------------------------------------------------------------------------------------|--|---------------------------------------------------------------------------|
| Sensitivity - $\mu V/V/mmHg$ 5.0 Nominal                                                                                                     |  | Temperature Coefficient of Sensitivity - $\%/^{\circ}C$ 0 $\pm$ 0.1       |
| Zero Balance - mmHg 0 $\pm$ 50                                                                                                               |  | Temperature Coefficient of Offset - mmHg/ $^{\circ}C$ 0 $\pm$ 0.3 Nominal |
| Zero Drift - mmHg < 2 in 4-hrs. after 5-min. warm-up                                                                                         |  | Operating Pressure - mmHg -30 to 300                                      |
| Excitation Voltage Range - Volts 4 - 8 (DC to 5 kHz)                                                                                         |  | Overpressure Withstand - mmHg -400 to 4000                                |
| Phase Shift < 5°                                                                                                                             |  | 15% Transducer Bandwidth - Hz > 200                                       |
| Excitation Impedance - Ohms 630 Nominal                                                                                                      |  | Symmetry 1 $\pm$ 5%                                                       |
| Signal Impedance - Ohms 300 $\pm$ 10%                                                                                                        |  | Defibrillator Withstand 5 in 5-min. at 360J                               |
| Calibration Requirement - mmHG 100 $\pm$ 3                                                                                                   |  | Risk Current of 115 VAC at 60 Hz - $\mu A$ < 5                            |
| Accuracy - mmHg $\pm$ 1 mmHg plus $\pm$ 1% of reading in -30 to 50 mmHg range, $\pm$ 1 mmHg plus $\pm$ 3% of reading in 50 to 300 mmHg range |  | Light Sensitivity - mmHg < 1 shift                                        |
| Compensated Temperature Range - $^{\circ}C$ 15 to 40                                                                                         |  | Mechanical Shock Withstand 3 falls from 1m                                |
| Storage Temperature Range - $^{\circ}C$ -25 to 70                                                                                            |  | Monitor Compatibility all monitors                                        |

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LogiCal, Sure-Cal, C-Fusor, Clear-Cuff, designmerket Medex og designmerket Smiths Medical er varemerker som tilhører Smiths Medical. Symbolet® angir at varemerket er registrert ved U.S. Patent and Trademark Office og i visse andre land. Alle andre navn og merker som er nevnt, er varemerker, handelsnavn og servicemerker som tilhører de respektive eiere.