

# Product Data Sheet

**Product Name: Cuff Pressure Gauge**

PEDIATRIC



SENSITIVE



POCKET



MONITOR



UNIVERSAL



## Order Information

REF	Description	Box
54-02-001	PEDIATRIC Cuff Pressure Gauge with hook	1
54-03-001	SENSITIVE Cuff Pressure Gauge with hook	1
54-04-000	POCKET Cuff Pressure Gauge with hook	1
54-05-000	MONITOR Cuff Pressure Gauge with hook	1
54-07-000	UNIVERSAL Cuff Pressure Gauge with hook	1

## Intended Use

REF 54-02-001:

The Cuff Pressure Gauge is used to adjust and control the cuff pressure of high-volume low-pressure cuffs of tracheal tubes.

REF 54-03-001/ REF 54-04-000/ 54-05-000/ REF 54-07-000:

The Cuff Pressure Gauge is used to adjust and control the cuff pressure of high-volume low-pressure cuffs of tracheal tubes and cuffs of supraglottic airways.

## Key data

REF	54-02-001	54-03-001	54-04-000	54-05-000	54-07-000
<b>Patient target group</b>	children	children and adults	children and adults	children and adults	children and adults
<b>Pressure range</b>	0-60 cmH <sub>2</sub> O	0-120 cmH <sub>2</sub> O	0-120 cmH <sub>2</sub> O	0-120 cmH <sub>2</sub> O	0-120 cmH <sub>2</sub> O
<b>Application</b>	Tracheal Tubes	Tracheal Tubes and Supraglottic Airways	Tracheal Tubes and Supraglottic Airways	Tracheal Tubes and Supraglottic Airways	Tracheal Tubes and Supraglottic Airways
<b>Box weight</b>	approx. 233 g	approx. 234 g	approx. 229 g	approx. 258 g	approx. 254 g
<b>Box dimensions (l) x (w) x (h) [mm]</b>	approx. 210 x 70 x 81	approx. 210 x 70 x 81	approx. 157 x 96 x 73	approx. 165 x 96 x 87	approx. 165 x 96 x 87
<b>Box content</b>	1 x Cuff Pressure Gauge 1 x Connecting Tubing for Cuff Pressure Gauge, length 100 cm 1 x Instruction for Use				

## Product Characteristics

<b>Sterile packed:</b>	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	<b>RoHs compliant:</b>	Yes <input type="checkbox"/>	No <input type="checkbox"/> N/A <input checked="" type="checkbox"/>
<b>Single use product:</b>	Yes <input checked="" type="checkbox"/> Connecting Tubing	No <input checked="" type="checkbox"/> Cuff Pressure Gauge	<b>REACH SVHC included:</b>	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>
<b>Phthalates</b> requiring labelling acc. to CLP Regulation (EC) 1272/2008	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>	<b>Use of natural rubber latex:</b>	Yes <input type="checkbox"/>	No <input checked="" type="checkbox"/>
<b>Shelf-life:</b>	5 years <input checked="" type="checkbox"/> Connecting Tubing	8 years <input checked="" type="checkbox"/> Cuff Pressure Gauge	<b>MRI suitable:</b>	Yes <input checked="" type="checkbox"/> Connecting Tubing	No <input checked="" type="checkbox"/> Cuff Pressure Gauge

Revision: 02/09.2021

**VBM Medizintechnik GmbH**

Einsteinstrasse 1 | 72172 Sulz a. N. | Germany  
Tel.: +49 7454 9596-0 | Fax: +49 7454 9596-33  
e-mail: info@vbm-medical.de | www.vbm-medical.de