

EU PRODUCT DATA SHEET CELOX GAUZE

Description:

Recommended by the US Committee on Tactical Combat Casualty Care (CoTCCC)¹ and is used by NATO forces and first responder units around the world. Celox Gauze is designed to stop bleeding from arterial injuries, road traffic accidents and other life-threatening injuries.

Presentation:

CE: 1.5m long x 7.6cm wide 'Z-fold' gauze

Features:

- Significantly reduces blood loss²
- Stops severe traumatic bleeding
- Works independently of the body's clotting mechanism
- Stops hypothermic bleeding and increases chance of survival³
- Celox gauze demonstrated a statistically significant improvement in survival (7%) and most evidently in the more severely injured⁴

How to use:

Tightly pack to source of the bleeding, above skin level. Apply a firm compression for three minutes, or until bleeding stops.

Indications for use:

To be used by trained emergency responders in the pre-hospital setting for temporary treatment of emergency life-threatening bleeding.

Patient Target Group: Adults and children, excluding neonates and infants.

Ordering Information:

FG Code: 5ft Z-Fold FG08838031 CE
100 packs per case

Storage:

Celox products should be stored in dry conditions at ambient temperature. No special storage conditions are required or indicated on the product labelling.

Shelf life five years from date of manufacture.



Packet Size: (Approx.)	CE
	Width - 11.1cm
	Height - 0.5cm
	Length - 13.2cm
Shipper Size: (Approx.)	Weight - 32g
	Width - 40cm
	Height - 16cm
	Length - 50cm
	Weight - 3.9kg

**FOR MORE INFORMATION
VISIT CELOXMEDICAL.COM**