



Su Instruments Pvt. Ltd	DEVICE	Thumb Tissue Forceps		Ref	l_1 ± 5.00	b_1 ± 0.5	b_2 ± 1.00						
	REF.NO.	SUBBD5577		SUBBD5577	145.00	3.25	8.75						
Capt. Shukar Deen Street Sialkot-51300-PAKISTAN	Material	Stainless Steel 410											
	Size	14.5cm											
Tel: +92-3006111353 Email: info@su-instruments.com www.su-instruments.com	Charge	MM/YY/FF		Page									
	Date	17-03-2025	Approved by	Director QC									

Material

Stainless Steel AISI 410

Legend

- 1: Space for labeling
- 2: Teeth
- 3: Surface profile of the manufacturer's choice

Applied Standards

DIN 58239-1: 2009-01

DIN 58298:2005-12 Medical instruments-materials, testing and testing

DIN 96042 Surgical Instruments Packing, Technical Specification

Application Areas

Tissue Forceps Thumb Tissue Forceps are used for grasping, holding, and manipulating various types of tissues in multiple procedures. ...

Others have greater tension and can hold thick, heavy tissue.

Thumb tissue forceps are commonly toothed in order to grip the tissue more firmly than serrated tips alone.

[illegible]