

Reference(s)	See Table Over Page
Description	Sterile, single use, rigid, transparent suction devices for use with suction connecting tubing. For general purpose surgical suctioning.
Intended Use	Pennine's Surgical Suction Probe range of devices are hand operated, single use, sterile devices. They are used as a conduit to remove bodily fluids and secretions, surgical tissue debris and irrigation fluids. They are used in surgical sites and body cavities to facilitate observation and/or to clear an obstruction. (GMDN 35917)
Duration of Use	Normally intended for continuous use of less than 60 minutes, but will be patient specific depending on the duration of the surgical procedure.
Product Information	Available with or without vacuum control (VC) vent Short mini, Mini fine, Midi and Maxi Yankauers available with smooth rounded end or Super Tip Short & Midi Yankauer available with Atrauma Tip to reduce tissue grab.
Materials	PVC (DEHP Free) Latex Free
Quality System	ISO 13485
EU Classification	Class IIa
Box Qty	See Table Over Page
Packaging	Sterile, double wrapped individual peel pouches
Sterilised by	Ethylene Oxide
Shelf Life	5 years



### Storage and Safety:

Toxicity	Manufactured from biocompatible materials.
Storage	Store in a clean, cool, dry place away from direct sunlight

© 2020 Pennine Healthcare or one of its subsidiaries.

All rights reserved. The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.

Product Code	Type	Length	CH	Qty. Per Box	UDI Code
LY-3400	Short Mini Tip	190mm	15ch	100	05033241019024
LY-3401	Short Mini Tip VC	190mm	15ch	100	05033241019055
LY-4400	Short Mini Supertip	195mm	15ch	100	05033241019642
LY-4401	Short Mini Supertip VC	195mm	15ch	100	05033241019666
LY-4420	Short Mini Super Tip (2 side	195mm	15ch	100	05033241019680
LY-3500	Mini Fine Tip	270mm	12ch	50	05033241019093
LY-3501	Mini Fine Tip VC	270mm	12ch	50	05033241019123
LY-4500	Mini Fine Supertip	275mm	12ch	50	05033241019703
LY-4500/NL	Mini Fine Super Tip 10° (No SIDE EYES)	275mm	12ch	50	05033241019772
LY-4501	Mini Fine Supertip VC	275mm	12ch	50	05033241019758
LY-4501/NL	Mini Fine Super Tip VC 0° (No	275mm	12ch	50	05033241019772
LY-3600	Midi Tip 30°	280mm	18ch	50	05033241019192
LY-3601	Midi Tip VC	280mm	18ch	50	05033241019246
LY-3600/PR	Midi Tip 10°	270mm	18ch	50	05033241019222
LY-4600	Midi Supertip	280mm	18ch	50	05033241019826
LY-4600/PR	Midi Super Tip 10°	280mm	18ch	50	05033241019857
LY-4601	Midi Supertip VC	280mm	18ch	50	05033241019888
LY-3700	Maxi Tip	275mm	22ch	50	05033241019390
LY-3700/F	Maxi Tip (with blue connector)	275mm	22ch	50	05033241019420
LY-3700/RF/FF	Maxi Tip (with straight connector)	275mm	22ch	50	05033241019451
LY-3701	Maxi Tip VC	275mm	22ch	50	05033241019482
LY-3702	Maxi Tip (No side eyes)	275mm	22ch	50	05033241019512
LY-4700	Maxi Supertip	275mm	22ch	50	05033241019925
LY-4700/PR	Maxi Super Tip 10°	275mm	22ch	50	05033241019949
LY-4701	Maxi Supertip VC	275mm	22ch	50	05033241019970



© 2020 Pennine Healthcare or one of its subsidiaries.

All rights reserved. The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.

Product Code	Type	Length	CH	Qty. Per Box	UDI Code
LY-3800	Maxi Plus Tip	280mm	28ch	50	05033241019529
LY-3801	Maxi Plus Tip VC	280mm	28ch	50	05033241019550
LY-3803	Maxi Plus Tip VC (No side eyes)	280mm	28ch	50	05033241019581
LY-3602	Atrauma tip, short soft	280mm	18ch	50	05033241019277
LY-3603	Atraumatic tip, short soft VC	280mm	18ch	50	05033241019291
LY-3606	Atrauma tip, short rigid	280mm	18ch	50	05033241019321
LY-3607	Atraumatic tip short rigid VC	280mm	18ch	50	05033241019338
LY-3608	Atraumatic tip, long rigid	280mm	18ch	50	05033241019352



Maxi Plus (four eyes)



Soft Rigid



Short Rigid



Long Rigid

© 2020 Pennine Healthcare or one of its subsidiaries.

All rights reserved. The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.

