

SPECIFICATION SHEET

MEETING YOUR RADIAL ACCESS NEEDS

Specially designed introducer kits, with or without Terumo hydrophilic polymer coating Mcoat™, for transradial access. The features of the Radifocus Introducer Kits II ensure smooth and easy access. The kink-resistant design helps maintain optimal vessel access and device passage throughout the entire procedure.

KEY FEATURES

SUTURE EYE

Aids in securing the sheath to skin

COLOUR-CODED SHEATH SUPPORT

For size determination

KINK-RESISTANT DESIGN

Helps keep the lumen open throughout the procedure

SMOOTH GUIDEWIRE-TO-DILATOR AND DILATOR-TO-SHEATH TRANSITIONS

For smooth penetration

MCOAT HYDROPHILIC COATING*
Designed for easy insertion and removal

SNAP-ON / CLICK-OFF DILATOR

Prevents backout during insertion

UNIQUE TERUMO CROSS-CUT VALVE

Avoids bleeding and air aspiration

RADIFOCUS™ *Introducer II Mcoat™*
Hydrophilic Coated Introducer Sheath

RADIFOCUS™ *Introducer II*
For radial access

* On RADIFOCUS™ Introducer II Mcoat™

ORDERING INFORMATION

Radifocus™ Introducer II – Mcoat™

Mcoat™ A-Kit

(sheath, dilator, 0.025" plastic mini guidewire, plastic IV catheter, syringe and scalpel)

Inner diameter	Sheath length		
	10 cm	16 cm	25 cm
5 Fr	RM*AF5J10SQW	RM*AF5J16SQW	RM*AF5J25SQW
6 Fr	RM*AF6J10SQW	RM*AF6J16SQW	RM*AF6J25SQW

Mcoat™ R-Kit

(sheath, dilator, spring mini guidewire and metallic entry needle)

Inner diameter	Guidewire compatibility	Sheath length		
		10 cm	16 cm	25 cm
5 Fr	0.021" (0.53 mm)	RM*RF5F10PQ	RM*RF5F16PQ	RM*RF5F25PQ
5 Fr	0.025" (0.64 mm)	RM*RF5J10PQ	RM*RF5J16PQ	RM*RF5J25PQ
6 Fr	0.021" (0.53 mm)	RM*RF6F10PQ	RM*RF6F16PQ	RM*RF6F25PQ
6 Fr	0.025" (0.64 mm)	RM*RF6J10PQ	RM*RF6J16PQ	RM*RF6J25PQ

Radifocus™ Introducer II – For radial access

(sheath, dilator, spring guide wire and metallic entry needle)

	Sheath length 7 cm		Sheath length 10 cm		
	Compatible with guidewire				
	Inner diameter	0.018" (0.46 mm)	0.025" (0.64 mm)	0.018" (0.46 mm)	0.021" (0.53 mm)
4 Fr	RS*R40A07PQ	RS*R40G07PQ	RS*R40A10PQ	RS*R40D10PQ	RS*R40G10PQ
5 Fr	RS*R50A07PQ	RS*R50G07PQ	RS*R50A10PQ	RS*R50D10PQ	RS*R50G10PQ
6 Fr	RS*R60A07PQ	RS*R60G07PQ	RS*R60A10PQ	RS*R60D10PQ	RS*R60G10PQ
7 Fr	-	-	-	RS*R70D10PQ	-



Terumo Corporation
+81 3 3374 8111

Terumo Europe NV
+32 16 38 12 11

Terumo Interventional Systems EMEA
+33 1 47 16 09 30

EMEA Sales Offices

Terumo Europe NV Africa Business Division
+32 16 38 13 08

Terumo Europe NV Benelux Sales Division
Belgium:
0800 14468
The Netherlands:
0800 0231938

Terumo Europe NV Emerging Market Division
+32 16 38 12 11

Terumo Deutschland GmbH
+49 6196 8023 0

Terumo Deutschland GmbH Österreich
+43 2236 378020

Terumo Deutschland GmbH Switzerland
+41 56 419 10 10

Terumo Europe España SL
+34 902 10 12 98

Terumo France S.A.S.
+33 1 30 96 13 00

Terumo Italia S.r.l.
+39 06 94 80 28 00

Terumo Russia LLC
+7 495 988 4740

Terumo Sweden AB
+46 3174 85 880

Terumo Middle East FZE
+971 4 292 0200

Terumo UK Ltd
+44 1276 480440

Terumo BCT Tıbbi Cihazlar Dağıtım ve Hizmetleri A.Ş.
+90 216 645 92 00

Terumo Poland Sp. z o.o.
+48 22 120 16 00

® Registered Trademark
Published by Terumo Europe NV