

EU PRODUCT DATA SHEET CELOX RAPID

Description:

The haemostatic of choice by the UK Ministry of Defense¹, Celox Rapid is our fastest-acting haemostatic gauze that removes the delay by speeding up packing time² and reducing compression time to 60 seconds, or until bleeding stops – typically saving over two minutes of treatment time, compared to other brands.

Presentation:

CE: 1.5m long x 7.6cm wide 'Z-fold' gauze.

Features:

- Can work with 60 seconds compression or until bleeding stops³
- Faster packing time than other leading brands²
- Rapid action reduces blood loss⁴
- Stops hypothermic bleeding⁵ as found in severe traumatic injuries

How to use:

Tightly pack to source of the bleeding, above skin level. If bleeding persists, apply a firm compression for 60 seconds or until bleeding stops.

Indications for use:

To be used by trained emergency responders in the pre-hospital setting for temporary treatment of emergency life-threatening bleeding.

Patient Target Group: Adults and children, excluding neonates and infants.


Ordering Information:

FG Code: FG08839021 CE
100 packs per case

Storage:

Celox products should be stored in dry conditions at ambient temperature. No special storage conditions are required or indicated on the product labelling.

Shelf life five years from date of manufacture.



Packet Size: (Approx.)	CE
	Width - 11.1cm
	Height - 0.5cm
	Length - 13.0cm
Shipper Size: (Approx.)	Weight - 30g
	Width - 40cm
	Height - 16cm
	Length - 50cm
	Weight - 3.9kg

**FOR MORE INFORMATION
VISIT CELOXMEDICAL.COM**