

# Product Data Sheet

## Product Name: Rebreathing Bag

| REF   |  |
|---|--|
| 66-00-011<br>66-00-015<br>66-00-012<br>66-00-013<br>66-00-014 |  <p>22 mm I.D.</p>                                 |
| 66-00-041<br>66-00-045<br>66-00-042<br>66-00-043<br>66-00-044 |  <p>22 mm I.D.</p> <p>Tube connection and cap</p> |
| 66-00-441<br>66-00-445<br>66-00-442<br>66-00-443<br>66-00-444 |  <p>22 mm I.D.</p> <p>5 mm hole</p>              |

## Order Information

| REF       | Description   | Box |
|-----------|---|-----|
| 66-00-011 | Rebreathing Bag, reusable, 500 ml standard, with 22mm I.D. cuff and loop                      | 1   |
| 66-00-015 | Rebreathing Bag, reusable, 1000 ml standard, with 22mm I.D. cuff and loop                     | 1   |
| 66-00-012 | Rebreathing Bag, reusable, 1500 ml standard, with 22mm I.D. cuff and loop                     | 1   |
| 66-00-013 | Rebreathing Bag, reusable, 2300 ml standard, with 22mm I.D. cuff and loop                     | 1   |
| 66-00-014 | Rebreathing Bag, reusable, 3000 ml standard, with 22mm I.D. cuff and loop                     | 1   |
| 66-00-041 | Rebreathing Bag, reusable, 500 ml with open tube connection and cap, 22mm I.D. cuff and loop  | 1   |
| 66-00-045 | Rebreathing Bag, reusable, 1000 ml with open tube connection and cap, 22mm I.D. cuff and loop | 1   |
| 66-00-042 | Rebreathing Bag, reusable, 1500 ml with open tube connection and cap, 22mm I.D. cuff and loop | 1   |
| 66-00-043 | Rebreathing Bag, reusable, 2300 ml with open tube connection and cap, 22mm I.D. cuff and loop | 1   |
| 66-00-044 | Rebreathing Bag, reusable, 3000 ml with open tube connection and cap, 22mm I.D. cuff and loop | 1   |
| 66-00-441 | Rebreathing Bag, reusable, 500 ml with 5 mm hole, 22mm I.D. cuff and loop                     | 1   |
| 66-00-445 | Rebreathing Bag, reusable, 1000 ml with 5 mm hole, 22mm I.D. cuff and loop                    | 1   |
| 66-00-442 | Rebreathing Bag, reusable, 1500 ml with 5 mm hole, 22mm I.D. cuff and loop                    | 1   |
| 66-00-443 | Rebreathing Bag, reusable, 2300 ml with 5 mm hole, 22mm I.D. cuff and loop                    | 1   |
| 66-00-444 | Rebreathing Bag, reusable, 3000 ml with 5 mm hole, 22mm I.D. cuff and loop                    | 1   |

## Intended Use

Rebreathing Bags are collapsible gas containers which are used in different anesthetic systems. They act as gas storage reservoir for the ventilator to balance discrepancy between patient's respiratory demand and gas flow delivered from the machine. The Rebreathing Bag also serves as a monitoring device for respiration during sedation. Furthermore, the Rebreathing Bag also permits manual ventilation and acts as a visual or tactile indicator of spontaneous breathing.

## Key data

|   |                              |                              |                              |                              |                              |
|---|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|
| <b>REF</b>  | 66-00-011                    | 66-00-015                    | 66-00-012                    | 66-00-013                    | 66-00-014                    |
| <b>Size</b>                                       | 500 ml                       | 1000 ml                      | 1500 ml                      | 2300 ml                      | 3000 ml                      |
| <b>Material</b>                                   | Silicone                     | Silicone                     | Silicone                     | Silicone                     | Silicone                     |
| <b>Colour</b><br>Bag<br>Cuff                      | translucent<br>green         | translucent<br>green         | translucent<br>green         | translucent<br>green         | translucent<br>green         |
| <b>Tube connection<br/>and cap</b>                | No                           | No                           | No                           | No                           | No                           |
| <b>5 mm hole</b>                                  | No                           | No                           | No                           | No                           | No                           |
| <b>Number of<br/>maximum<br/>reprocessings</b>    | 100                          | 100                          | 100                          | 100                          | 100                          |
| <b>Box weight</b>                                 | 95 g                         | 107 g                        | 118 g                        | 175 g                        | 160 g                        |
| <b>Box dimension</b><br>Length<br>Width<br>Height | 21.5 cm<br>15.0 cm<br>3.5 cm |

|   |                              |                              |                              |                              |                              |
|---|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|
| <b>REF</b>  | 66-00-041                    | 66-00-045                    | 66-00-042                    | 66-00-043                    | 66-00-044                    |
| <b>Size</b>                                       | 500 ml                       | 1000 ml                      | 1500 ml                      | 2300 ml                      | 3000 ml                      |
| <b>Material</b>                                   | Silicone                     | Silicone                     | Silicone                     | Silicone                     | Silicone                     |
| <b>Colour</b><br>Bag<br>Cuff                      | translucent<br>green         | translucent<br>green         | translucent<br>green         | translucent<br>green         | translucent<br>green         |
| <b>Tube connection<br/>and cap</b>                | Yes                          | Yes                          | Yes                          | Yes                          | Yes                          |
| <b>5 mm hole</b>                                  | No                           | No                           | No                           | No                           | No                           |
| <b>Number of<br/>maximum<br/>reprocessings</b>    | 100                          | 100                          | 100                          | 100                          | 100                          |
| <b>Box weight</b>                                 | 97 g                         | 110 g                        | 120 g                        | 177 g                        | 165 g                        |
| <b>Box dimension</b><br>Length<br>Width<br>Height | 21.5 cm<br>15.0 cm<br>3.5 cm |

|   |                              |                              |                              |                              |                              |
|---|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|
| <b>REF</b>  | 66-00-441                    | 66-00-445                    | 66-00-442                    | 66-00-443                    | 66-00-444                    |
| <b>Size</b>                                       | 500 ml                       | 1000 ml                      | 1500 ml                      | 2300 ml                      | 3000 ml                      |
| <b>Material</b>                                   | Silicone                     | Silicone                     | Silicone                     | Silicone                     | Silicone                     |
| <b>Colour</b><br>Bag<br>Cuff                      | translucent<br>green         | translucent<br>green         | translucent<br>green         | translucent<br>green         | translucent<br>green         |
| <b>Tube connection<br/>and cap</b>                | No                           | No                           | No                           | No                           | No                           |
| <b>5 mm hole</b>                                  | Yes                          | Yes                          | Yes                          | Yes                          | Yes                          |
| <b>Number of<br/>maximum<br/>reprocessings</b>    | 100                          | 100                          | 100                          | 100                          | 100                          |
| <b>Box weight</b>                                 | 95 g                         | 107 g                        | 118 g                        | 175 g                        | 160 g                        |
| <b>Box dimension</b><br>Length<br>Width<br>Height | 21.5 cm<br>15.0 cm<br>3.5 cm |

## Product Characteristics

|   |                              |  |
|---|------------------------------|--|
| <b>Sterile packed:</b>  | Yes <input type="checkbox"/> | No <input checked="" type="checkbox"/> |
| Sterilization method:   | n/a                          |  |
| <hr/>   |                              |  |
| <b>Single use product:</b>  | Yes <input type="checkbox"/> | No <input checked="" type="checkbox"/> |
| <hr/>   |                              |  |
| <b>Phthalates</b> requiring<br>labelling acc. to CLP<br>Regulation (EC) 1272/2008 | Yes <input type="checkbox"/> | No <input checked="" type="checkbox"/> |
| <hr/>   |                              |  |
| <b>Shelf-life:</b>  | 5 years                      |  |

|   |   |  |
|---|---|--|
| <b>RoHS compliant:</b>                  | Yes <input type="checkbox"/>            | No <input type="checkbox"/>            |
|   | n/a <input checked="" type="checkbox"/> |  |
| <hr/>                                   |   |  |
| <b>REACH SVHC included:</b>             | Yes <input type="checkbox"/>            | No <input checked="" type="checkbox"/> |
| <hr/>                                   |   |  |
| <b>Use of natural rubber<br/>latex:</b> | Yes <input type="checkbox"/>            | No <input checked="" type="checkbox"/> |

---

Revision: 02/03.2021

---

**VBM** Medizintechnik GmbH

Einsteinstrasse 1 | 72172 Sulz a. N. | Germany  
Tel.: +49 7454 9596-0 | Fax: +49 7454 9596-33  
e-mail: info@vbm-medical.de | www.vbm-medical.de