



Medline Industries Ltd
3rd Floor, Quayside Wilderspool Business
Park
Greenalls Avenue
Warrington WA4 6HL
United Kingdom

Tel.: +44 844 334 5237
Fax: +44 844 334 5238
www.medline.eu/uk



Technical Sheet

Secretion bag

Reference(s)

References	Volume [ml]	Tubing length [cm]	Feature	Sterility	Pack Factor
690515	1500	90	without NRV	Non-sterile	200
690520	2000	90	without NRV	Non-sterile	200
690415	1500	90	with NRV	Non-sterile	200
690420	2000	90	with NRV	Non-sterile	200

Description

Secretion Bag for collecting discharged body fluids, such as urine from foley catheters. The Bag is facilitated with an approximately graduation to offer visibility on urine output. Not intended for precise measurement. The Bag offers fluid collection of up to 1,5l or 2l depending on the version. It is equipped with a gradual step connector and comes in 2 different versions, with or without non-return valve (NRV). The tube comes in a length of 90cm and an outer diameter of 6.5mm (+/- 0.15mm). As the product is not equipped with a drainage valve, the secretion bag needs to be replaced after filling. Indwelling time up to 7 days. Not to be used in cases of in case of known hypersensitivity to the materials used

Main materials

Component	Material
Bag	PVC
Tube	PVC (MT-202)
Step connector	PVC
Protective Cap	PE

Pack factor

Refer to the reference table

Packaging

polybag non-sterile

Legal Manufacturer

CONOD Medical Co., Limited
No. 11 Hongfeng Road, Baimao Industrial Park,
Guli Town 215532 Changshu City, Jiangsu Province
PEOPLE'S REPUBLIC OF CHINA
In compliance with EN ISO 13485

Regulation

In compliance with Medical Devices Regulation 2017/745
MD Class: I non-sterile

•Please refer to complete labeling before using the product.

Other Requirements

Shelf Life: 5 years
Storage: +5°C to +30°C, in a dry place, away from sunlight

© 2023 Medline Industries LP or one of its subsidiaries. All rights reserved
TDS_Secretion bag non-sterile_690xxx_EN01

Retention period: End of life and support period



Latex Free